




## Helps to limit the inhalation of microorganisms and their transmission by air

### Specifications:

Product: small-sized particle filtering mask  
Model name: C-PEF99  
Class: FFP2 NR  
Expanded size: 12x11 cm  
Executive standard:  
EN 149:2001 +A1 2009  
Shelf life: 3 years from the date of manufacture  
Storage: 41°F / 70°F (5°C / 20°C), humidity 10 % / 85%



### Usage:

Protects children's nose and mouth from inhaling infectious agents.



Product code: **FFP2x2KD**  
Paraf code: **X85L5K 981649585**  
EAN code: **7640173131267**

The best protection  
for children